Viztek technical publications Document Number OPAL -2006-DC3

OPAL-RAD

DICOM Conformance Statement

TABLE OF CONTENTS

TABLE OF CONTENTS	2
A.0 Introduction	3
A.1 Implementation Model	3
A. 1.1 Application Data Flow Diagram	4
A.1.2 Functional Definition of Application Entities	5
A. 1.3 Sequencing of Real World Activity	6
A.2 Application Entity specifications	7
A.2.1 Application Entity specifications for OPAL -RAD Server	7
A. 2.1.1 Verification	7
A. 2.1.2 Storage	7
A. 2.1.3 Storage Commitment	8
A. 2.1.4 Query / Retrieve	8
A. 2.2 Application Entity specifications for OPAL-RAD Acquire	10
A. 2.2.1 Verification	10
A. 2.2.2 Storage	10
A.2.3 Application Entity specifications for OPAL-RAD Film Acquire	11
A. 2.3.1 Verification	11
A 2.3.2 Storage	11
A.2.4 Application Entity specifications for OPAL-RAD Professional Workstation	11
A. 2.4.1 Storage	11
A. 2.4.2 Query / Retrieve	11
A. 2.4.3 Printing	12
A. 2.4.4 Media Storage	12
A. 2.5 DICOM Conformance for RIS Gateway Interface	12
A. 2.5.1 Detached Patient Management	12
A. 2.5.2 Detached Study Management	13
A. 2.5.3 Detached Study Component Management	15
A. 2.5.4 Detached Results Management	15
A. 2.5.6 Detached Interpretation Management	16
A. 2.5.7 Structured Reporting	17
A.3 Association Establishment Policies	19
A. 3.1 General	19
A. 3.2 Number of Associations	19
A. 3.3 Asynchronous Nature	19
A. 3.4 Implementation Identifying Information	19
A.3.5 Association Initiation / Acceptance Policy	19
A. 3.6 Associated Real World Activity	19
A. 3.7 Proposed Presentation Context Table	19
A.4 Communication Profiles	20
A. 4.1 Supported Communication Stacks	20
A. 4.1.1 <i>OSI Sta</i> ck	20
A. 4.1.2 TCP/ IP Stack	20
A. 4.1.3 Point to Point Stack	20
A.5 Extensions / Specializations / Privatizations	20
A.6 Configuration	20
A.6.1 AE Title / Presentation Address Mapping	20
A.7 Support of Extended Character Sets	20

A.O INTRODUCTION This document provides DICOM conformance information for the OPAL -R AD Telemedicine PACS produced by Viztek, Incorporated, and assumes the reader is familiar with the components and terminology intrinsic to DIC OM 3.0 protocol. The Service Classes and Information Objects supported by OPAL-R AD are described here in full detail, in accordance with Part 2 of the DIC OM standard. OPAL-R AD uses DICOM services to import, transfer and store DICOM images for analysis and processing.

A.1 IMPLEMENTATION MODEL OPAL-R AD exchanges data between the Service Class User (SCU) and Service Class Provider (SCP) over any TC P/IP network. OPAL -R AD Server (SC P) can communicate with three SC U components - OPAL-RAD Professional Workstation, OPAL-R AD Film Acquire and OPAL-RAD Acquire.

A.1.1 Application Data Flow Diagram Application data flow revolves around one SCP and three SCU components in OPAL-R AD. The data flow diagram for a fully implemented OPAL -R AD PACS is shown in Figure A.1. A functional OPAL -RAD system must include the Server (SC P) component and at least one SC U component (Professional Workstation, Film Acquire or Acquire). Interactions between components and outside sources are described. Only DIC OM functioning is listed inside each OPAL-R AD component process. The OPAL -R AD Professional Workstation is implemented with two separate processes — ePatient and Diagnostic Viewer. Professional Workstation can initiate and receive query/retrieve associations, while Server can only receive them.




Figure A.1 Application Data Flow Diagram

A.1.2 Functional Definition of Application Entities

- OPAL-RAD Server (SCP) is executed when the machine is powered on. It waits for an association request from a remote SCU client and initiates associations in the following ways.
- 1. QUERY: Waits for a C-FIND-RQ from a remote AE. The Server will send a C-FIND-R SP message with the matching information (or none if the query is unsuccessful). The remote AE c n cancel the query by sending a C-CANCEL-FIND-RQ message to OPAL-RAD Server.
- 2. RETRIEVE: Waits for a C-MOVE -RQ or C-GET -RQ from remote AE. The Server responds to the requesting AE using either C-MOVE-RSP or C-G ET-RSP (respective to the request) and sends the requested information on a separate association.
- 3. VERIFICATION: Waits for an association from remote AE to perform echo tests. If C-ECHO-RQ is received, a C-ECHO-RSP is send back with the status message "success."
- 4. STORAGE: A DICOM storage commitment is initiated if Server receives a C-STORE-R Q from OPAL-R AD Acquire, Film Acquire, or attached DICOM compliant modality. If the request is accepted, a patient study and series is stored in an Oracle database and a C-STORE -RSP is sent to the perspective AE.
- OPAL-RAD Acquire (SC U) is executed at the discretion of the user.
- 1 .STORAGE: Acquire initiates association with OPAL-RAD Server and sends a storage request (C-STORE RQ) to save a DICOM study and series. If the association is made and request is accepted, the data is transferred to the Server for storage.
- OPAL -RAD Film Acquire (SCU) is executed at the discretion of the user.
- 1 .STORAGE: Film Acquire initiates association with OPAL-RAD Server and sends a storage request (C-ST ORE-R Q) to save a DICOM study and series. If the association is made and request is accepted, the data is transferred to the Server for storage.
- OPAL –RAD Professional Workstation (SCU) is executed at the discretion of the user.
- 1. QUERY: Workstation initiates an association with a remote AE to query images on a remote SCP or SCU. An initial Patient-Lev el query is sent with a C-FIND-RQ operation to obtain a patient list based on a set of search criteria. Afterward, the Workstation can initiate a subsequent Study-Lev el C-FIND-RQ query for that particular patient. After performing a Study-Level query, the Workstation can initiate a Series-Level query for that particular patient and study (issuing another C-FIND-RQ). Professional Workstation can cancel a query by sending a C-CA NCEL- FIND-RQ message to the remote AE.
- 2. PRINT: Initiates association with DIC OM compliant print devices. If association is accepted, Professional Workstation sends a DICOM Print request (C-PRINT -RQ) and respective data to the device.
- 3. VERIFICATION: Initiates an association with a remote AE for performing an Echo request. When a successful association is made, it sends a C-E CHO-RQ on the same association. If a remote AE makes a successful association with Workstation and a C-ECHO-RQ is receive d, the Workstation sends a C-ECHO-RSP to the remote AE along with the message "success."
- 4. RETRIEVE: Initiates an association with two remote AEs (Source and Destination). After a successful Query to the Source AE node, a C-MOVE-R Q is sent to the source AE, supplying the destination AE Title. The Source AE then replies with a C-MOVE-RSP, transferring the requested information from the Source AE to the Destination AE. 5. MEDIA STORAGE: Can initiate extraction of DICOM information from a DICOM DIR structure on any local or removable media that is connected to the Workstation.

A.1.3 Sequencing of Real World Activity

There are several real world activities that initiate DICOM associations with OPAL-RAD Server (SCP) or other AEs. The number of DIC OM associations held between OPAL -RAD Server and remote AEs is only limited by the system's resources. Real world activities are listed below.

A.1.3.1 OPAL-RAD Server

No real world actions are triggered directly from OPA L -RAD Server or OPAL-R AD Server's interface. Real world actions are initiated only through OPAL-RAD Professional Workstation, Acquire, Film Acquire, a DICOM compliant modality, or another remote

A. A.1.3.2 OPAL-RAD Acquire

Storage – An association is initiated with OPAL-RAD Server (SCP) when the user presses "Save" from the Acquire main interface. A DICOM stud y is then created and sent to the SCP, which stores the DIC OM in a database.

A.1.3.3 OPAL-RAD Film Acquire Storage – An association is initiated with OPAL-RAD Server (SCP) when the user presses "Save" from the Film Acquire main interface. A DICOM stud y is then created and sent to the SCP, which stores the DIC OM in a database.

A.1.3.4 OPAL-RAD Professional Workstation

Viewer initiates associations from several real world activities.

- 1. Query When the user clicks "Update Patients" in the Workstation's DICOM Console screen, a DICOM association is created with the respective SCP or SC U. An association is also created when the user presses "Update Studies" or "Update Series" for each respective list in the DICOM Console Screen.
- 2. Retrieve When the user has loaded a Source and Destination AE node in the DICOM Console screen, a move operation can be conducted by pressing the "Move" button. The user must perform a "Query" operation with the Source AE prior to this. If the request is successful, the selected series, study or patient is copied and transferred from the source AE to the destination AE.
- 3. Print An association is made with any DICOM compliant printer when the user selects "DICOM Print" from the Print menu option in the Diagnostic Viewer screen.
- 4. Echo If the user clicks on "Echo Src" (meaning Echo Source) or "Echo Dest" (meaning Echo Destination) from the DICOM Console screen, a "Verification" association is made with the respective AE node. The user can also initiate an association from selecting "Echo" from the "Print" menu option in the Diagnostic Viewer screen. The user must enter an Application Entity and press the button labeled "Echo" to make the association.
- 5. Media Storage Professional Workstation can extract information from a DICOMDIR structure by clicking on the "Import DICOM DIR" button located on the standard toolbar in the Diagnostic Viewer Screen. This can also be initiated by selecting "Import DICOMDIR" from the "File" menu option. The user is then prompted to enter a drive/storage location or "browse" local and removable directories for the desired structure.

A.2 APPLICATION ENTITY SPECIFICATIONS

A.2.1 Application Entity specifications for OPAL-RAD Server

OPAL-RAD Server supports the following Application Entities in the SCP role.

A.2.1.1 Verification

so	P Class Name	SOP Class UID
Veri	ificatio n	1.2.840.10008.1.1

A.2.1.2 Storage

SOP Class Name	SOP Class UID
CR Im age Storage	1.2.840.10008.5.1.4.1.1.1
Digi X ray Pres Image Storage	1.2.840.10008.5.1.4.1.1.1.1
Digi X ray Proc Image Storage	1.2.840.10008.5.1.4.1.1.1.1.1
Digi Mammo Pres Image Storage	1.2.840.10008.5.1.4.1.1.1.2
Digi Mammo Proc Image Storage	1.2.840.10008.5.1.4.1.1.1.2.1
Digi Intra Oral Pres Image Storage	1.2.840.10008.5.1.4.1.1.1.3
Digi Intra Oral Proc Image Storage	1.2.840.10008.5.1.4.1.1.3.1
CT Im age Storage	1.2.840.10008.5.1.4.1.1.2
US M ulti Image Storage (Ret)	1.2.840.10008.5.1.4.1.1.3
MR Im age Storage	1.2.840.10008.5.1.4.1.1.4
NM Im age Storage(Ret)	1.2.840.10008.5.1.4.1.1.5
US I mage Storage (Ret)	1.2.840.10008.5.1.4.1.1.6
NM Im age Storage	1.2.840.10008.5.1.4.1.1.20
US Im age Storage	1.2.840.10008.5.1.4.1.1.6.1
US M ulti Image Storage	1.2.840.10008.5.1.4.1.1.3.1
SC Im age Storage	1.2.840.10008.5.1.4.1.1.7
RT Im age Storage	1.2.840.10008.5.1.4.1.1.481.1
RT Do se Storage	1.2.840.10008.5.1.4.1.1.481.2
RT Structure Set Storage	1.2.840.10008.5.1.4.1.1.481.3
RT Treatm ent Record Storage	1.2.840.10008.5.1.4.1.1.481.4
RT Plan Storage	1.2.840.10008.5.1.4.1.1.481.5
X-Ray Angiographic Image Storage	1.2.840.10008.5.1.4.1.1.12.1
X-R ay Radiofluoroscopic Image Storage	1.2.840.10008.5.1.4.1.1.12.2
X-Ray Angiographic Bi-Plane Image Storage	1.2.840.10008.5.1.4.1.1.12.3

	,
StandabneOverlay Storage	1.2.840.10008.5.1.4.1.1.8
Standabne Curve Storage	1.2.840.10008.5.1.4.1.1.9
Waveform WV Mod ality Storage	1.2.840.10008.5.1.4.1.1.9.1
Waveform ECG Mod ality Storage	1.2.840.10008.5.1.4.1.1.9.1.1
Waveform AU Mo dality Storage	1.2.840.10008.5.1.4.1.1.9.1.2
Waveform Hemodynamic Storage	1.2.840.10008.5.1.4.1.1.9.2.1
Standalone Modality LUT St orage	1.2.840.10008.5.1.4.1.1.10
StandabneVOILUT Storage	1.2.840.10008.5.1.4.1.1.11
Gray Softcopy Pres StateStorage	1.2.840.10008.5.1.4.1.1.11.1
VLE ndoscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.1
VL Microsc opic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.2
VL Sli de Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.3
V L Photographic I mage Storage	1.2.840.10008.5.1.4.1.1.77.1.4
SR Basic Text	1.2.840.10008.5.1.4.1.1.88.11
SR Enhanced	1.2.840.10008.5.1.4.1.1.88.22
SR Com prehensive	1.2.840.10008.5.1.4.1.1.88.33
PET Im age Storage	1.2.840.10008.5.1.4.1.1.128
PET C urve Storage	1.2.840.10008.5.1.4.1.1.129
-	

A.2.1.3 Storage Commitment

SOP Class Na me	SOP Class UID
Storage Commitment Push Model Class	1.2.840.10008.1.20.1
Storage Commitment Push Model Instance	1.2.840.10008.1.20.1.1
Storage Commitment Pull Model Class	1.2.840.10008.1.20.2
Storage Commitment Pull Model In stance	1.2.840.10008.1.20.2.1

A.2.1.4 Query/Retrieve

SOP Class Na me	SOP Class UID
Patient Root QR Find	1.2.840.10008.5.1.4.1.2.1.1
Patient Root QR Move	1.2.840.10008.5.1.4.1.2.1.2
Patient Root QR Get	1.2.840.10008.5.1.4.1.2.1.3
Study Root QR Find	1.2.840.10008.5.1.4.1.2.2.1
Study Root QR Move	1.2.840.10008.5.1.4.1.2.2.2
Study Root QR Get	1.2.840.10008.5.1.4.1.2.2.3
Patient/StudyRootQR Find	1.2.840.10008.5.1.4.1.2.3.1

Patient/Study Root QR Move	1.2.840.10008.5.1.4.1.2.3.2
Patient/Study Root QR Get	1.2.840.10008.5.1.4.1.2.3.3

A.2.2 Application Entity specifications for OPAL-RAD Acquire

OPAL -RAD Acquire supports the following Application Entities in the SC U role.

A.2.2.1 Verification

SOP Class Name	SOP Class UID
Verificatio n	1.2.840.10008.1.1

A.2.2.2 Storage

SOP Class Name	SOP Class UID
SC Im age Storage	1.2.840.10008.5.1.4.1.1.7

A.2.3 Application Entity specifications for OPAL-RAD Film Ac quire OPAL-RAD Film Acquire supports the following Application Entities in the SCU role.

A.2.3.1 Verification

SOP Class Name	SOP Class UID
Verificatio n	1.2.840.10008.1.1

A.2.3.2 Storage

SOP Class Name	SOP Class UID
SC Im age Storage	1.2.840.10008.5.1.4.1.1.7

A.2.4 Application Entity specifications for OPAL-R AD Professional Workstation OPAL -RAD Professional Workstation supports the following Application Entities in the SCU role.

A.2.4.1 Storage

SOP Class Name	SOP Class UID
Verificatio n	1.2.840.10008.1.1

A.2.4.2 Quer y / Retrieve

SOP Class Na me	SOP Class UID
Patient Root QR Find	1.2.840.10008.5.1.4.1.2.1.1
Patient Root QR Move	1.2.840.10008.5.1.4.1.2.1.2
Patient Root QR Get	1.2.840.10008.5.1.4.1.2.1.3
Study Root QR Find	1.2.840.10008.5.1.4.1.2.2.1
Study Root QR Move	1.2.840.10008.5.1.4.1.2.2.2
Study Root QR Get	1.2.840.10008.5.1.4.1.2.2.3
Patient/Study Root QR Find	1.2.840.10008.5.1.4.1.2.3.1
Patient/Study Ro ot QR Move	1.2.840.10008.5.1.4.1.2.3.2
Patient/Study Root QR Get	1.2.840.10008.5.1.4.1.2.3.3

A.2.4.3 Printing

SOP Class Na me	SOP Class UID
Basic Film Session	1.2.840.10008.5.1.1.1
Basic Fil m Bo x	1.2.840.10008.5.1.1.2
Basic Gray Image Box	1.2.840.10008.5.1.1.4
Basic Color Im age Box	1.2.840.10008.5.1.1.4.1
Basic Annotation Box	1.2.840.10008.5.1.1.15
Reference Image Box	1.2.840.10008.5.1.1.4.2
Print Job	1.2.840.10008.5.1.1.14
Printer	1.2.840.10008.5.1.1.16
Printer Instance	1.2.840.10008.5.1.1.17
VOIL UT B ox	1.2.840.10008.5.1.1.22
Presentation LUT	1.2.840.10008.5.1.1.23
Image Overlay Box	1.2.840.10008.5.1.1.24
Print QueueInstance	1.2.840.10008.5.1.1.25
Print Queue Management	1.2.840.10008.5.1.1.26
Stored Print Storage	1.2.840.10008.5.1.1.27
Hardcopy Gray Image Storage	1.2.840.10008.5.1.1.29
Hardcopy Color Image Storage	1.2.840.10008.5.1.1.30
Basic Gray Print Management Meta	1.2.840.10008.5.1.1.9
Basic Co lor Print Management Meta	1.2.840.10008.5.1.1.18
Reference Gray Print Management Meta	1.2.840.10008.5.1.1.9.1
Reference Color Print Management Meta	1.2.840.10008.5.1.1.18.1

A.2.4.4 Media Storage

SOP Class Name	SOP Class UID
Media Storage Dir Storage	1.2.840.10008.1.3.10

A.2.5 DICOM Conformance for RIS Gateway Interface OPAL-RAD DICOM RIS Gateway supports the following Application Entities in the SCU and SCP role. All N-E V ENT and SCU/SCP roles are supported for all service classes per the DI COM 2000 standard.

A.2.5.1 Detached Patient Management

Attribute Name	Tag
Specific CharacterSet	(0008,0005)
Patient Name	(0010,0010)
Patient ID	(0010,0020)
Issuer of Patient ID	(0010,0021)
Other Patient Ids	(0010,1000)
Other Patient Names	(0010,1001)
Patient Telephone Numbers	(0010,2154)
Patient Address	(0010,1040)
Patient Bir th Date	(0010,0030)
Patient Sex	(0010,0040)
Patient W eight	(0010,1030)
Ethnic Group	(0010,2160)
Patient Religious Preference	(0010,21F0)
Patient Data Confidentiality Constraint Desc.	(0040,3001)
Patient State	(0038,0500)
Pregnancy Status	(0010,21C0)
Medical A lerts	(0010,2000)
Contrast Allergies	(0010,2110)
Special Needs	(0038,0050)
Referenced Study Sequence	(0008,1110)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)
Referenced Visit Sequence	(0008,1125)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)
Reference Patient Alias Sequence	(0038,0004)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)

A.2.5.2 Detached Study Management

Attribute Name	Tag
Specific CharacterSet	(0008,0005)
Study ID	(0020,0010)
Study ID Issuer	(0032,0012)
AccessionNumber	(0008.0050)
Study Instance UID	(0020,000D)
Study Status ID	(0032,000A)
Study Priority ID	(0032,000C)
Scheduled Study Start Date	(0032,1000)
Scheduled Study Start Time	(0032,1001)
Scheduled Study Stop Date	(0032,1010)
Scheduled Study Stop Time	(0032,1011)
Scheduled Study Lo cation	(0032,1020)
Scheduled Study Lo cation Application Entity Title	(0032,1021)
Requesting Service	(0032,1033)
Requesting Physician	(0032,1032)
Requested Procedure Description	(0032,1060)
Requested Procedure Code Sequence	(0032,1064)
>Code V alue	(0008,0100)
>Coding Scheme Designator	(0008,0102)
>Code Meaning	(0008,0104)
Study A rrival Date	(0032,1040)
Study A rrival Ti me	(0032,1041)
Study Date	(0008,0020)
Study Time	(0008,0030)
Study Completed Date	(0032,1050)
Study Completed Time	(0032,1051)
Study Verified Date	(0032,0032)
Study VerifiedTime	(0032,0033)
Study Read Date	(0032,0034)
Study Read Time	(0032,0035)
Name of Physician(s) Reading Study	(0008,1060)
Reason For Study	(0032,1030)

Referenced Patient Sequence	(0008,1120)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)
Referenced Visit Sequence	(0008,1125)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)
Referenced Results Sequence	(0008,1100)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)
Referenced Study Component Sequence	(0008,1111)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)

A.2.5.3 Detached Study Component Management

Attribute Name	Tag
Specific CharacterSet	(0008, 0005)
Referenced Study Sequence	(0008,1110)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UI D	(0008,1155)
Retrieve Application Entity Title	(0008,0054)
Modality	(0008,0060)
Study Description	(0008,1030)
A cquistion In Study	(0020,1004)
Study Status ID	(0032,000a)

A.2.5.4 Detached Results Management

Attribute Name	Tag
Specific CharacterSet	(0008,0005)
Results ID	(4008,0040)

Results ID Issuer	(4008,0042)
Impressions	(4008,0300)
Referenced Study Sequence	(0008,1110)
> Referenced SO P Class UID	(0008,1150)
> Referenced SOP In stance UTD	(0008,1155)
Referenced Interpretation Sequence	(4008,0005)
> Referenced SO P Class UID	(0008,1150)
> Referenced SOP In stance UTD	(0008,1155)

A.2.5.6 Detached Interpretation Management

Attribute Name	Tag
Specific C haracter Set	(0008,0005)
Interpretation ID	(4008,0200)
Interpretation ID Issuer	(4008,0202)
Interpretation Type ID	(4008,0210)
Interpretation Status ID	(4008,0212)
Interpretation Recorded Date	(4008,0100)
Interpretation Recorded Time	(4008,0101)
Interpretation Recorder	(4008,0102)
Interpretation Transcription Date	(4008,0108)
Interpretation Transcription Time	(4008,0109)
Interpretation Transcriber	(4008,010A)
Interpretation Author	(4008,010C)
Interpretation Tex t	(4008,010B)
Referenced Results Sequence	(0008,1100)
> Referenced SOP Class UID	(0008,1150)
> Referenced SOP In stance UID	(0008,1155)
Interpretation Approver Sequence	(4008,0111)
>Interpretation Approval Date	(4008,0112)
>Interpretation Approval Time	(4008,0113)
>Physicians Approving Interpretation	(4008,0114)

A.2.5.7 Structured Reporting

SR Document Series Module Attributes

Attribute Name	Tag
Modality	(0008,0060)
Series Instance UID	(0020,000E)
Series Number	(0020,0011)
Referenced Study Component Sequence	(0008,1111)
>Referenced SOP Class UID	(0008,1150)
> Referenced SOP In stance UID	(0008,1155)

SR Document General Module Attributes

Attribute Name	Tag
Instance Number	(0020,0013)
Completion Flag	(0040,A4 91)
Completion Flag Description	(0040,A4 92)
VerificationFlag	(0040,A4 93)
Content Date	(0008,0023)
Content Time	(0008,0033)
Verifying Observer Sequence	(0040,A073)
>Verifying Observer Name	(0040,A075)
>V erifying Observer Identification C ode Sequence	(0040,A088)
>V erifying Organization	(0040,A027)
>Verification DateTime	(0040,A030)
Predecessor Documents Sequence	(0040, A 360)
Identical Documents Sequence	(0040,A525)
Referenced Request Sequence	(0040,A370)
>Study Instance UID	(0020,000D)
>Referenced Study Sequence	(0008,1110)

>>Referenced SOP Class UID	(0008,1150)
>>Referenced SOP Instance UID	(0008,1155)
>Accession Number	(0008,0050)
>Placer Order Number/Imaging Service Request	(0040,2016)
>Filler Order Number/Imaging Service Request	(0040,2017)
>Requested Procedure ID	(0040,1001)
>Requested Procedure Description	(0032,1060)
>Requested Procedure Code Sequence	(0032,1064)
Performed Procedure Code Sequence	(0040,A372)
Current Requested Procedure Evidence Sequence	(0040,A375)
Pertinent Other Evidence Sequence	(0040,A385)

A.3 ASSOCIATION ESTABLISHMENT POLICIES

- A.3.1 General There is no limit on maximum PDU size. It is usually limited by the available system memory. The default value is 16 K (1 6,384 byte).
- A.3.2 Number of Associations There is no inherent limitation in the number of Associations accepted by OPAL RAD Server, Professional Workstation, Acquire and Film Acquire. System resources usually limit it. Server, Professional Workstation, Acquire and Film Acquire initiate only one Association at a time.
- A.3.3 Asynchronous Nature Asynchronous operations are not supported. A request for asynchronous transfer is rejected.
- A.3.4 Implementation Identifying Information The Implementation Class UI D is "1.2.840.114062" and Implementation Version Name is "OPAL-R AD v4".
- A.3.5 Association Initiation/Acceptance Policy Both OPAL-RAD Acquire and professional Workstation initiate a new association only with operator interaction. For instance, the operator may send a study to a film printer from a Workstation. The SCP does not initiate a new association. It only processes the association from a remote node, optionally alerting the operator about the activity. It places no limitations on who may connect to it for verification purposes. The remote AE must exist and be configured properly in the Server's database to perform query and retrieve operations with OPAL-R AD.
- A.3.6 Associated Real World Activity All real world activity for initiating an association is conducted through the SCU applications, Professional Workstation, Acquire and Film Acquire. Real World Activities are described in A.1.3.

A.3.7 Proposed Presentation Context Table OPAL -RAD accepts the following transfer syntaxes.

Name List	UID List
JPEG Baseline	1.2.840.10008.1.2.4.50
JPEG Exte nded	1.2.840.10008.1.2.4.51
JPEG Lossless No Pred	1.2.840.10008.1.2.4.57
JPEG Lossless	1.2.840.10008.1.2.4.70
RLE Lossless	1.2.840.10008.1.2.5
Implicit VR Little En dian	1.2.840.10008.1.2
Explicit VR Little En dian	1.2.840.10008.1.2.1
Explicit VR Big Endian	1.2.840.10008.1.2.2

If a remote client proposes multiple transfer syntaxes to OPAL-R AD Server, then preference is given to the transfer syntax that matches the host machine's architecture. The maximum number of Transfer Syntaxes presented in an Association is 64.

A.4 COMMUNICATION PROFILES

A.4.1 Supported Communication Stacks OPAL-RAD provides TCP/IP Communication Support as defined in Part 8.

A.4.1.1 *OS I Sta*ck OSI stack is not supported.

A.4.1.2 TCP/IP Stack TC P/IP stack is inherited from the operating system, and is supported by OPAL-R AD.

A.4.1.3 Point to Point Stack Point to Point stack is not supported.

A.5 EXTENSIONS/SPECIALIZATIONS/PRIVATIZATIONS None.

A.6 CONFIGURATION - All configurable parameters are set through OPAL-R AD Server at runtime. The configuration is password protected at Administrator level.

A.6.1 AE Title / Presentation Address Mapping Local AE title for Acquire is "OPAL Acquire SCU". Local AE title for Film Acquire is "OPAL Film Ac SCU". Local AE title for Professional Workstation is "OPAL Viewer SC U". Local AE title for Server is "OPAL Server SC U". Node IP address, net mask, hostname, hostname aliases are configured by the local system administrator. TCP/IP port is configurable. Default is 104 DICOM.

A.7 SUPPORT OF EXTENDED CHARACTER SETS No extended character sets are supported.