


VIZTEK

U-ARM BASIC



Your Practice. Your Preference.
Your Productivity Solution.

VIZTEK U-ARM BASIC

The Viztek U-Arm Basic offers outstanding image quality with easy positioning. The stand is fully motorized with dual speeds and also has automatic table and upright positioning.

The Basic has multiple 17x17 panel options available and multiple generator sizes available. The system works on the highly advanced, Viztek developed, Ultra software that allows for the friendliest user experience.



VIZTEK U-ARM BASIC



Pediatric Imaging Initiative (optional)

Viztek's Pediatric Imaging software package is designed to assist in identifying dose reduction opportunities.

The user-defined settings enables a reminder based on age to remove the grid to reduce dose. Patient age is determined in the software by the date of birth and prompts x-ray technique variation based on suggested dose per age group.



2217 US Highway 70 East • Garner, NC 27529
Phone: 800.366.5343
Email: sales@viztek.net • www.viztek.net

Ultra Software Features

- User-friendly controls
- Foldering allows for greater control and organization
- Improved workflow
- Single screen access to editing
- Windows 8/7 capability



- Table Weight - 96lbs
- Table Capacity - 400lbs

Optional Features

- Automatic Exposure Control (AEC)
- Weight Bearing Stand (WBS)
- Dose Area Product meters (DAP)

