


# Digital Radiography for the Mobile World.



## ViZion™ + DR Mobile

The ViZion™ + DR mobile kit features the lightweight wireless panel along with the reliable SR130 portable X-ray system. The digital package uses a mounted laptop tray and convenient panel slide to hold the laptop and the panel. The wireless panel is available in both Gadolinium (GOS) and Cesium (CSI) technologies.

### ViZion™ + Wireless DR Specifications

Weight: 7.7 lbs (GOS) 7.9 lbs (CSI)  
Carbon Fiber and Mg Alloy Construction  
Pixel Matrix - 2304 x 2800  
3.4 lp/MM (GOS) 3.7 lp/mm (CSI)  
Cassette-sized  
Lithium Ion Battery (2 Included)

### SR 130 X-ray System Specifications

3.0kw Output - Fully Regulated mA / kVp  
KVP Range 40-100  
MA Range 15 and 30  
MAS Range 0.15-120

ViZion is a trademark of Konica Minolta Healthcare Americas, Inc.  
© 2017 Konica Minolta Healthcare Americas, Inc.