


WIRELESS
VIZION DR +



Your Practice. Your Preference.
Your Productivity Solution.

VIZION DR + WIRELESS CSI

Viztek's ViZion DR + wireless panel provides exceptional image quality paired with advanced software features and functionality. The 14x17 panel uses enhanced Cesium technology and has unique features like drop notification emails.

It is ideal for use in both new systems as well as retrofit rooms. The versatile ViZion DR + wireless panel utilizes the ViZion Acquire software for image capture, manipulation, and sending. Viztek's acquire software offers useful tools to enhance the image acquisition process.

Panel Specifications

Weight - 8.1lbs
Pixel Size - 150um
Pixel Array - 2288x2784
3.4 lp/mm
Dimensions - 384x460x15.2mm
A/D Conversion 14 Bits

ViZion Software Features

Generator Integration
Patient Review
Reverse Black/White
Multi-Scale Contrast Enhancement
and Dynamic Range
Optimized Image Processing
Parameters



2217 US HWY 70 East • Garner, NC 27529 • Phone: 800.366.5343
Email: info@viztek.net • www.viztek.net

Viztek is a trademarks of Viztek LLC. ©2016 Viztek LLC. All rights reserved.

SA-151A Rev A