



## Viztek Opal CR System

Viztek Opal-CR is the fast and easy way to convert your facility to digital today! The easy-to-use tabletop CR can produce 56 plates per hour while offering crisp image quality. Reading images at 10 pixels/mm, the compact CR can be used in any environment, including mobile units. Viztek Opal CR is the fastest way to take existing analog imaging into the digital realm. Additionally, the CR also offers significant speed and image quality over analog film processing.

### System Features

- Single Cassette Feed
- 2 Cassettes Included
- Opal Front-End Viewing Software
- Data Acquisition : 20 bits/pixel
- Pixel Matrix Size : 3480x4248
- 14 x 17 cassette size : 10 pixels/mm
- Power Requirement : 120V/60Hz
- Dimension : 27.2 x 27.6 x 18.2"



2217 US HWY 70 EAST • Garner, NC 27529  
Phone: 800.366.5343 • Fax: 919.772.2530  
Email: [info@viztek.net](mailto:info@viztek.net) • [www.viztek.net](http://www.viztek.net)

## CR BENEFITS

- Available in Multiple Configurations
  - Available with Opal Mini PACS Bundle
  - Additional Viewers Available
- Opal-Acquisition Software
- Market Leading Image Quality
- Compact Table-top Convenience
- Supports DICOM Modality Worklist
- Plates available in four sizes
- Allows for Importing Outside CDs
- Burn Patient CDs with Viewer
- Scan Documents in to Exams\*
- Fast Installation and Applications Process
- Training and 24/7 support for hardware and software components
- Eliminate the cost associated with processing film
- Mobile Kits Available

\* Requires Document Scanner - Not Included