

LEGGERA DR SOLUTIONS

Leggera DR Solutions are available in a cesium wireless cassette size wifi panel and a 4143 fixed panel. Both offer superior image quality, low dose, and offer images within seconds of exposure for speedy throughput.

The 35x43cm wireless panel features an ultra-light design and is highly versatile, fitting in various x-ray room environments. At only 6.17lbs and with up-to 8-hours of battery life, the (ISO 4090) cassette size panel uses ethernet communication technology and includes a back-up connection cable. An optional battery charger is available and can charge up to three batteries at a time.

The fixed 4143 panel is an ideal solution for various x-ray room configurations including the FIT Kit. This incorporates the Fixed and Wireless panel in a dual room.

| Leggera DR Solutions | FIXED | WIRELESS |
|-------------------------|------------------|------------------------|
| Pixel Size | 148um | 148um |
| Dimensions | 19.7x19.3x1.8 in | ISO 4090 Cassette Size |
| Weight | 27.55lbs | 6.17lbs |
| DQE @ 0lp/mm | 64% | 66% |
| AD Conversion | 16 Bits | 16 Bits |
| Image Display Time | Less than 5 sec | Preview in 3 sec |
| X-Ray Gen Voltage Range | 40-150kv | 40-150kv |



