

Your Practice. Your Preference Your Productivity Solution.

straight arm or

Viztek's Straight Arm DR offers an advanced digital x-ray solution in an efficient and compact design. The system features a 17 x 17 detector for superior image quality. It allows for all imaging views that are required in radiology including wheelchair and table work for maximum patient flexibility. The Straight Arm DR maximizes every square foot of space to avoid overcrowded rooms and provides more space for staff to maneuver.

Specifications

- Swivel arm range of rotation: 135 degrees with detents on 30,0,90,120 degrees
- Swivel arm vertical movement: 39 inches
- Minimum floor to detector center distance: 27.6'
- Detector Angulations: 45/-45 degrees
- Manual Collimator
- DR adaptation for 17x17 and 14x17 application
- Available with CSI or GOS panels
- Different Generator Sizes Available
- Weight Bearing Stand is Optional



2217 US HWY 70 East • Garner, NC 27529 • Phone: 800.366.5343 Email: info@viztek.net • www.viztek.net raight arm of